

USE MANUAL OF MAGIC CHAIR

Magic Chair Manual



When operating this product, please read the manual of this instrument carefully

Content

1.Instrument introduction	3
Principle&advantages.....	3
2.Precautions	3
2.1 Operation/safety precautions	4
3.Instrument installation.	6
3.1 component	6
4.Operation process of the instrument.....	7
4.1.Starting interface	7
5.Maintenance and repair	9
6.Common troubleshooting method.....	9
7.Responsibility of the manufacturer.....	10
8.Storage and transportation management.....	10
9.Warranty terms	11
10.Technical Parameters.....	11
11.Parts List.....	12

1.Introduction

1.1 Instrument principle

The "Magic Chair" subverts the traditional treatment of urinary incontinence. Through the revolutionary technology "HIFEM" high-intensity focused electromagnetic field, it stimulates motor neurons to control the pelvic floor muscles, and triggers repeated hyper-contraction motions of the pelvic floor muscles to rebuild the strength and endurance of the pelvic floor muscles and the urinary system.

1.2 Instrument advantages

"Magic Chair" can activate the pelvic floor muscle tissue and promote collagen regeneration; effectively enhance the pelvic floor muscles, restore firmness and improve elasticity; strengthen the pelvic floor muscles, fight off incontinence, urine leakage, dripping, lower abdomen protrusion, lower pelvic relaxation and other problems; get rid of difficult orgasm and other problems. The pelvic floor muscles are triggered to tighten approximately 12,000 times within 30 minutes, which effectively strengthens the pelvic floor muscles, improves urine penetration, and relaxes in private areas. The effect far exceeds the general Kegel exercise.

- Using the most advanced non-invasive vaginal shrinking technology, no need to undress, safe, comfortable, and no abstinence period.**
- A course of treatment is better than completing 10,000 Kegel exercises, and the treatment effect is more significant.**
- Each treatment only needs to sit for 30 minutes, with zero pain, zero anesthesia, and zero recovery period.**

2. Precautions

After the instrument is opened, check whether the accessories are complete according to the instrument configuration list. If you have any problems, please contact the seller.

2.1 Operation taboos

- 1、 It's forbidden for pregnant women, breastfeeding women, and menstruating women.
- 2、 Those who are close to menstruation should be cautious, which may cause early menstruation or abdominal cramps.
- 3、 It's forbidden for those position which implant metal or electronic equipment (pacemaker, defibrillator, neurostimulator, drug pump, etc.).
- 4、 It's forbidden for those with major diseases such as malignant tumors and cancer.
- 5、 It's forbidden for those with heart disease, arrhythmia, epilepsy, severe hypertension, and diabetes.
- 6、 It's forbidden for those who have artificial filling in the body.
- 7、 It's forbidden for those with major surgery patients, and it can also be used according to the recovery situation and after consulting the doctor.
- 8、 If there is a wound on the skin, avoid the wound or don't do it.
- 9、 Before use, the metal objects worn on the body must be removed. When using, avoid the metal objects on the clothes.
- 10、 For those people with too thick fat layer are not recommended to use this device. This device is mainly for muscle tissue training. For those people with too thick fat layer, the energy hardly reach the muscle layer to cause muscle contraction and affect the effect. It 's recommended to do it after lipolysis.

2.2 Safety precautions

- 1、 When the power plug is dirty, wipe it clean in time, and make sure that there are no foreign objects before connecting to the socket.
- 2、 It must be connected high-power supply socket with three plug at 220-230v,50Hz, and remove the power cord when moving the equipment.
- 3、 Don't touch the power plug and pull the power cord after plugging in.
- 4、 Do not touch the power cord and plug with wet hands.
- 5、 Do not force the instrument or the parts around the power supply.
- 6、 Do not use the instrument at high temperatures, and do not place heating equipment near the instrument.
- 7、 After the instrument is started, do not open the power socket of the instrument at will, which will cause failure.

- 8、 If there is smoke or odor when using the instrument, please unplug the power supply and consult after-sales service.
- 9、 Do not touch the lips and eyes when the handle is activated to avoid injury.

2.3 Operation precautions

- 1、 Please read the manual carefully, and then use it according to the method in the manual. When the instrument is not used for a long time, please unplug the power supply.
- 2、 The instrument must be used after consulting with experts or after professional training.
- 3、 When using this instrument, first turn down the intensity, and then gradually adjust to the applicable intensity.
- 4、 After using it each time, you need to stand by for ten minutes before turning off the power. If you need to change the operator, please use it after 5 minutes.
- 5、 If there is any abnormality on the patient when using the instrument, please stop the operation immediately and consult an expert.
- 6、 If the instrument is abnormal, please confirm the manual again before operation; if the abnormality is not removed, do not disassemble the instrument by yourself, and contact after-sales service for consultation.

3.Installation


3.1 component




4. Operation process

4.1 Starting interface


Choose a language



Male / Female
Mode selection



Note: 1. Choose according to different genders and parts



2. The back, neck, and privacy can be independently selected, and the handle strength/time can be adjusted according to different customers, and then click the start/stop button to operate.

Note: 1 select the part that needs to be operated:

2 adjust the intensity/frequency time:

5. Maintenance and repair

5.1 Maintenance

5.1.1 Daily inspection

Managers must regularly inspect and control instruments and components.

Whether the power socket is reliably connected to the ground.

Before use, check whether the power cord of it is damaged, and stop using the instrument if it is damaged.

5.1.2 cleaning

Periodically cleaning the equipment to ensure the sanitation and good operation of it. You can wipe the outside of the device with a damp cloth when it's powered off, using soapy water or a neutral detergent.

5.1.3 Fuse replacement

When replacing the fuse, in order to prevent electric shock, you must turn off the power and unplug the power plug.

5.2 Repair

After purchasing the machine, you need to keep the warranty card and send it to the manufacturer or agent to have the warranty effect. From the date of purchase (by receipt), the company's product quality is guaranteed for one year free of charge. (Warranty with the purchase receipt and warranty card.) After the warranty period has expired, please contact the manufacturer directly to provide professional after-sales service. To repair this machine, you must contact the company's business department or dealer, and do not repair it without authorization. If you open the package without consent, the company may not provide warranty.

6. Common troubleshooting method.

Failure phenomenon	Causes	Exclusion and resolution
Equipment not working	The power plug is not inserted into the socket The power switch is not on Failure of power fuse	check power check the power switch replace the fuse

Alarm prompt code

When the treatment device malfunctions, it automatically stops the treatment work and sends a fault alarm.

Display fault codes. As shown in Figure 5.

Fault code list:

Code	Explanation	Notes
E3:	Abnormal coil temperature	
E1:	Abnormal temperature of the function box	

7. Responsibility of the manufacturer

- The products we produced should meet the product standards. If there is any non-compliance, the company is responsible for replacement and return.
- In case of personal and property damage caused by product manufacturing defects, our company shall bear corresponding liabilities, but the company shall not be liable for the following situations;
- The purchaser and operating user do not use in accordance with the operating specifications of the manual, and do not use in accordance with the safety precautions and warning instructions.
- Violating product transportation, installation, use, maintenance, storage, etc., to cause product, personal, or property damage.
- The company implements life-long maintenance and technical services for the products. If the manufacturing quality problems occur within the specified warranty period, they will be repaired for users free of charge
- When accessories need to be replaced, original accessories must be purchased from our company, otherwise the adverse consequences will be borne by the user.

8. Storage and transportation management

Storage: The instrument should be stored in the original packaging box, placed in a well-ventilated room, the packaging box should be high, the ambient temperature is $-10^{\circ}\text{C}\sim 40^{\circ}\text{C}$, the relative humidity is not more than 93%, and it is not allowed to be placed in harmful gas. Surrounding flammable, explosive materials and corrosive gases.

Transportation: The instrument in the packaged condition is suitable for road, railway, aviation and waterway transportation. During loading and unloading and transportation, it should be protected from severe vibration and shock, and should not be exposed to moisture, and cannot be mixed with flammable and corrosive substances. The specific requirements are as specified in the order contract.

9.Warranty terms

The warranty period of this instrument is 1 year (Note: except for artificial damage)

The warranty period of accessories is 3 months (Note: except for artificial damage)

If the instruments beyond the service life are to be used continuously, they shall be tested and approved by the engineer, and the performance of it shall be tested regularly by the engineer in the subsequent use. The scrapped instruments and accessories exceeding the service life shall be disposed according to the relevant local regulations.

Maintenance record sheet			
Product name		Purchase date	
Dealer name		Telephone	
Dealership address			
Customer name		Telephone	
Customer address			
Repair date	Failure	Repair content	

10.Technical Parameters

Item	Parameter
Power	110-220 V50 /60Hz
power consumption	2.3KVA
Electromagnetic wave energy	0~7T
Pulse	300μs (microsecond)
Case size (cm)	86*80*71(cm)
Packing Weight(kg)	63kg

11.Parts List

Serial No.	part	quantity	unit
1	Host (qualified product)	1	PCS
2	power cord	1	PCS
3	Warranty card	1	PCS
4	use manual	1	PCS